tousimis[®] 931 Touch Screen Display Overview


4.1 Start Up




NOTE:

Pressure measurement is not indicative of the amount of LCO₂ remaining in tank. The same "PSI" reading may be noted even after the liquid CO₂ drops below the level of the Syphon (Dip) Tube. Weight is the best way to properly indicate remaining LCO₂.


1. All 5 metering valves (COOL, SLOW FILL, FILL, BLEED, PURGE/VENT) have been factory preset (See Chapter 5 Illustration for Metering Valves). There is no need to adjust these Metering Valves. Open the main LCO₂ Tank valve. The Syphon (Dip) Tube LCO₂ Tank should have between 25 lbs - 50 lbs (12.5 kg -25 kg) net weight of LCO₂ in order to operate the 931. (Subtract the tare weight of the tank from the tank weight to get the net weight of LCO₂).
2. Turn power switch ON.



3. The 931 will start up with the following screen display.



4. The 931 is ready to operate when the following screen is displayed.



Press “**Enter User/Sample Information**” to enter User/Sample Information.
You may also press “**Main Menu**” to start operation without entering User & Sample Information (See #5 & #6).

If PROCESS CHAMBER is under pressure (above 50 PSI), the 931 will display “CHAMBER UNDER PRESSURE.”

Press VENT VALVE Switch to open VENT VALVE in order to start a new process run.

If the previous process run was interrupted by any unexpected reason such as a power outage during the process run, DO NOT OPEN VENT VALVE.
You may resume the rest of process run steps using “Manual Mode” from Main Menu screen (#7).

5. Enter User Name or ID (maximum 18 spaces long) and press **Enter**.

The screen displays the prompt "Please Enter User Name or ID" at the top. Below the prompt is a text input field. Underneath the input field is a numeric keypad with buttons for digits 1-0 and a "Clr" button. Below the numeric keypad is a QWERTY keyboard layout. At the bottom of the screen are three buttons: a blue "SKIP" button, a "Space" button, and an "Enter" button.

6. From the following screen, Enter Sample ID or Description (maximum 18 spaces long) and press **Enter**.

The screen displays the prompt "Please Enter Sample Description" at the top. Below the prompt is a text input field. Underneath the input field is a numeric keypad with buttons for digits 1-0 and a "Clr" button. Below the numeric keypad is a QWERTY keyboard layout. At the bottom of the screen are three buttons: a blue "SKIP" button, a "Space" button, and an "Enter" button.

7. Main Menu will then be displayed. Press any button to initiate desired mode.

The Main Menu screen features a dark header bar with the "tousimis" logo on the left, the text "MAIN MENU" in the center, and user/sample/date information on the right. Below the header, the current chamber status is shown as "0 PSI" and "12°C". The main area contains six large orange buttons with icons and labels: "AUTO" (checkmark icon), "MANUAL" (hand icon), "STASIS" (circular arrows icon), "RECIPE" (fork and knife icon), "NEW USER" (person icon), and "SETTINGS" (wrench icon). At the bottom, a status bar displays "CTRL: V.2.39.01", "GUI: v.0.9.304", and "© 2012 tousimis".


Annotations:

- Title/Current Mode:** Points to the "MAIN MENU" text in the header.
- User/Sample /Date Information:** Points to the text "Hanz Wetzler Tp 334 Btch 16653 01/26/2012 10:43 AM" in the header.
- Chamber Temperature/ Pressure Display Area:** Points to the "0 PSI" and "12°C" display.
- Control and Status Display Area:** Points to the six main control buttons.
- Software Version:** Points to the "CTRL: V.2.39.01" and "GUI: v.0.9.304" text in the footer.



4.2 Auto Mode

1. Factory Default Time and Pressure Values may be used or updated from the Auto Mode Setup Screen.




Press **NEXT** to advance screen or to modify preset values press appropriate buttons. For example, to modify “Slow FILL” time, press **Slow FILL** button; to update BLEED/VENT Pressure, press **BLEED/VENT** button.


Press **NEXT** button to start “Auto Mode.”

To restore tousimis® recommended preset values, press **FACTORY DEFAULT** button. Exit the mode by pressing **MAIN MENU** button.


2. To update Time, use arrow buttons or slide bar to change value and press **SAVE**. **UNDO** button restores tousimis® recommended default values for typical process run. Press **MAIN MENU** or **AUTO SETUP** button to exit without saving.



3. To update BLEED/VENT Pressure, use arrow buttons or slide bar to change value and press **SAVE**. **UNDO** button restores tousimis® recommended value for typical process run. Press **MAIN MENU** or **AUTO SETUP** button to exit without saving.



4. Press **NEXT** button from **AUTO MODE SETUP** page once Time/Pressure modifications are input to satisfaction.



NOTE:

* Use Ultrapure Alcohols only!
(i.e. I.P.A., Methanol, Ethanol).

Use minimum 99.5%+ purity Alcohols for best results.


NOTE:

DO NOT introduce Acetone into the Process Chamber unless Acetone Kit (#8770-91) has been previously installed.


NOTE:

Never use channel locks or pliers for tightening the knurl nuts. Always hand tighten!


5. Pour sufficient alcohol into process chamber (Typically 10-15ml) to ensure complete sample coverage. Next, transfer the sample holder into the open process-chamber. Make the transfer quickly in order to avoid the samples' exposure to air and moisture. The chamber O-Ring should lie completely in its groove. NEVER grease the O-Ring. Place the chamber lid down over the chamber using the 3 knurl nuts. Tighten each knurl nut hand tight with equal pressure. Once lid is secure press **COOL** button.




6. As the chamber temperature slowly begins to drop, you may hear the LCO₂ circulating through the unit. The 931 will continue cooling by itself until the chamber temperature reaches 0°C (±5°C).



7. At this cut off point, the cooling will automatically stop and display the following message:



8. Press the **SLOW FILL** button to advance. The 931 requires confirmation of chamber security to further continue process. Double check chamber security and press **CONFIRM**.







NOTE

PAUSE button will pause operation and give option to change preset value during process run. See Chapter 4.5 "Pause & Resume during Process Run."

9. The 931 will begin to fill the chamber with LCO₂. *From this point forward, the 931 will automatically advance through all process modes sequentially until completion.* During the SLOW FILL mode, the LCO₂ will enter and fill the Process Chamber slowly.




10. The 931 will advance to FILL Mode. The chamber temperature is automatically maintained below 10°C (±5°C). This is normal.



NOTE

Pressure OFF button may be used in case of poor chamber sealing. Chamber O-Ring is able to be replaced without sample exposure to the air.

* Pressing the **PRESSURE OFF** button will slowly reduce chamber pressure while allowing alcohol to remain in the Process Chamber (if necessary).






NOTE

The waste alcohol will exit the 931 Chamber via the PURGE/VENT Exhaust into the Carboy or Solvent Drain.

11. After the FILL mode expires, the 931 will automatically advance into the PURGE mode. The 931 will remain in the PURGE mode for the duration of process time.




12. Upon Completion of the PURGE mode, the unit will automatically advance into a POST-PURGE FILL mode in which the chamber fills with LCO₂.



NOTE

Consult Check-Out Data Sheet in the Appendix of this 931 User Manual.

13. Upon completion of the POST-PURGE FILL mode, the HEAT mode will activate. The HEAT mode is when the samples are carried through the "Critical Point." Both the pressure and temperature will steadily increase.







NOTE:

- + : CP Temperature
- : CP Pressure
- ⊕ : Critical Point


14. When the chamber pressure reaches and goes beyond 1072 PSI, it will stabilize in the range of 1350PSI ($\pm 5\%$ @ 20°C). As the temperature reaches 31-35°C, the unit has achieved the “critical point” and this is where the ‘tousimis equilibrium’* cycle starts. The + sign between the Pressure Temperature display indicates the Process Chamber has reached critical point temperature. The ○ sign indicates the Process Chamber has reached critical point pressure.



15. If chamber pressure has not reached critical pressure prior to critical temperature “CP Failure” will be displayed. Refer “Critical Point Failure Correction (Chapter 4.7)” section.



16. At the end of the ‘tousimis equilibrium’ period, the 931 will automatically advance into the BLEED mode.






NOTE:

Step #17 is optional and NOT necessary with each process run. This is more of a preventive maintenance detail.


17. At this point, you can measure the BLEED rate (if desired) using the 20 SCFH Flow Meter supplied, by attaching the Flow Meter to the outlet of the COOL/BLEED EXHAUST. The BLEED rate has been factory preset to decompress the chamber at a rate of 100-150 PSI/min.

The flow rate should read 9-12 SCFH at the onset of the BLEED mode. This setting should yield an average decompression rate of approximately 100-150 PSI/min reduction in pressure. This pressure reduction flow rate is the desired decompression rate between 1300 PSI to 550 PSI.

Once the BLEED rate has been adjusted, remove the Flow Meter and reconnect COOL/BLEED Exhaust Hose.




18. At the preset pressure (factory default: 550 PSI), the 931 will automatically advance from the BLEED mode into the VENT mode.



19. It is not necessary to re-adjust the PURGE-VENT metering valve flow rate as it is preset. The Process Chamber should then come to atmospheric pressure after approximately 3 minutes in this VENT mode.

20. When chamber pressure reaches less than 10 PSI, the 931 will display “Process Complete.”



21. At this point, the chamber lid may be removed by alternatively and evenly loosening the Knurl Nuts (Never attempt to 'force' open).

22. The sample(s) can then be removed from the Process Chamber. Reseal the Process Chamber with the chamber lid to help keep it clean and moisture free.

23. The next process run may be initiated provided there is sufficient LCO₂. Otherwise, the system power may be turned OFF.



4.3 Manual Mode

1. The “Manual Mode” enables the 931 to operate without a preset process run sequence. User may press buttons to advance each mode.
2. Prior to initiating the COOL mode, pour sufficient alcohol (Typically 10-15ml) into Process Chamber to cover your sample. Quickly transfer sample holder(s) into the Process Chamber.
The chamber O-Ring should lie completely in its groove. NEVER grease the O-Ring. Place the chamber lid down over the chamber using the 3 knurl nuts. Tighten each knurl nut hand tight with equal pressure. Once lid is secure press **COOL** button.

Actual Manual Mode Initial Screen Display




Secure Chamber




NOTE:

Never use channel locks or pliers for tightening the knurl nuts. Always hand tighten!




3. The 931 will continue cooling and maintaining the chamber temperature near 0°C during manual COOL Mode until next mode has been selected by user.



NOTE

Chamber Security is important to maintain chamber pressure. Chamber Lid must be properly sealed as per p.30 "Star Pattern" Chamber Closure Sequence.

4. Press the **SLOW FILL** button to advance the process. The 931 requires Confirmation of "Chamber Security" to continue the process. Double check "Chamber Security" and press **CONFIRM**.




NOTE

One may observe the chamber through the Chamber Lid Window to make sure LCO₂ covers sample completely.

5. The 931 will begin to fill the chamber with LCO₂. During the SLOW FILL Mode, the LCO₂ will enter and slowly fill the Process Chamber.




6. Press the **FILL** button to advance to FILL Mode.



NOTE

Pressure OFF button may be used in case of poor chamber sealing. Chamber O-Ring is able to be replaced without sample exposure to the air.


* Pressing the **PRESSURE OFF** button will slowly reduce chamber pressure, while allowing alcohol to remain in the Process Chamber (if necessary).




NOTE

PAUSE button will pause operation and give option to change preset value during process run, "Resume" operation or terminate the process run. See section 4.5 "Pause & Resume During Process Run."

7. Press **PURGE** button to advance into the PURGE Mode. During PURGE Mode the waste alcohol will exit the 931 Chamber into the Carboy or Solvent Drain.




8. The 931 will shut down after 15 minutes of PURGE Mode. Checking exhaust tubing at this point is recommended. The 931 can resume any operational mode by pressing the appropriate button.




NOTE

A two minute POST-PURGE FILL time is recommended or until the chamber liquid level is full.

9. Upon completion of the PURGE Mode, you may press the **FILL** button to accomplish POST-PURGE FILL Mode, during which the chamber only fills with LCO₂.



10. Press **HEAT** button to advance into the HEAT Mode. The HEAT Mode is the stage in which the samples are carried through the "Critical Point." Both the pressure and temperature will steadily rise as the Process Chamber has active HEAT.






NOTE:

- + : CP Temperature
- : CP Pressure
- ⊕ : Critical Point


11. When the Process Chamber pressure reaches and goes beyond 1072 PSI, it will stabilize in the range of 1350 PSI ($\pm 5\%$ @ 20°C). As the temperature achieves 31°C, the Process Chamber has achieved the “critical point”. The + sign on the Pressure/Temperature display indicates the Process Chamber has reached the critical point temperature, the ○ sign indicates the Process Chamber has reached the critical point pressure.




NOTE:

If needed, measure BLEED rate using "Auto Mode" BLEED Instructions (Chapter 4.2 Step #17).

12. It is recommended to wait for at least 2 minutes during the ‘tousimis equilibrium’ period prior to pressing **BLEED** button to advance into the BLEED mode. It is recommended to maintain BLEED Mode down to 550 PSI.



13. Press **VENT** button to advance VENT Mode to complete a process run.





4.4 Stasis Mode

(See Chapter 4 prior to proceeding)

1. Stasis Mode processes samples using the identical Process Steps as the “Auto Mode” except for an additional “Stasis Step” (see Chapter 4 Stasis Diagram).

The Stasis Time and the number of Stasis Cycle need to be setup (see #3) to operate “Stasis Mode”.

2. You may either use the Factory Default settings or Custom Setup the desired Time and Pressure Values for each Process Mode and press **NEXT** (See Chapter 4.2).



NOTE

With the “Stasis Mode”, the user is able to program the times that the sample will sit in STASIS. After the **PURGE** Mode the Post **PURGE** STASIS will repeat.



NOTE

Maximum Time for Slow FILL, FILL, PURGE and Post-PURGE Fill Mode is 15 minutes.
Maximum CP Time is 180 Minutes.

Navigation Buttons

Default Time/Pressure Values for each Process Mode

3. Using the Up/Down arrows, you may setup the desired “Stasis Times”, “Stasis Cycles” and “Stasis PURGE Time”. Press **START** to initiate the Stasis Mode.

Up/Down Arrows

Press **RESET** button to restore the factory default values.

Press either **MAIN MENU** or **AUTO SETUP** buttons to exit without saving values.

For the STASIS Process, Autosamdri®-931 will use two different PURGE Time. Normal PURGE Time will be used for the initial PURGE Process before the 931 advances into Stasis Cycle.

Stasis PURGE Time will be used during Stasis Cycle.

Most Alcohol will be exit from process chamber during Normal PURGE Time.

The Stasis Purge Time recommended is to set it to less than 30% of normal PURGE Time to help reduce LCO₂ consumption.

Example: STASIS Process Settings

Slow FILL: 3 Min

FILL: 6 Min

Normal PURGE Time: 10 Min

Post-PURGE FILL: 4 Min

Critical Point: 10 Min

Stasis Time: 2 Hours

Cycle: 3 Times

Stasis PURGE Time: 5 Min

Process Sequence Example:

COOL → SLOW FILL (3 MINS) → FILL (6 MINS) → PURGE TIME (10 MINS)
→ POST-PURGE-FILL (4 MINS)
→ STASIS #1 (2 HOURS) → COOL → SLOW FILL (3 MINS)
→ STASIS PURGE TIME (3 MINS) → POST-PURGE FILL (4 MINS)
→ STASIS #2 (2 HOURS) → COOL → SLOW FILL (3 MINS)
→ STASIS PURGE TIME (3 MINS) → POST-PURGE FILL (4 MINS)
→ STASIS #3 (2 HOURS) → COOL → SLOW FILL (3 MINS)
→ STASIS PURGE TIME (3 MINS) → POST-PURGE FILL (4 MINS)
→ HEAT → CRITICAL POINT (tousimis®-EQUILIBRIUM 10 MINS)
→ BLEED → VENT

NOTE:


Total PURGE Time = 19 Min

10 Mins (Normal PURGE Time) + 3 x 3 Mins (Stasis PURGE Time) = 19 Mins


4. Once the **START** button is pressed, the 931 will operate as per Auto Mode from this point forward until POST-PURGE FILL Mode.



Stand By Mode



PURGE Mode




POST-PURGE FILL Mode



NOTE

All valves leading into and out of the Process Chamber are closed to maintain LCO₂ inside of chamber during "STASIS Mode".

5. Upon completion of the POST-PURGE FILL Mode, the 931 will automatically advance into the first cycle of the STASIS Mode.




The sample will now stay submerged in LCO₂ allowing time for the intermediary fluid to dissipate from sample interior out into the surrounding LCO₂.

6. After the initial STASIS Mode is completed, the 931 will advance into COOL Mode and resume the process.





NOTE

It is normal to see the chamber temperature to gradually increase during Stasis Mode.



7. The 931 will then automatically advance into the SLOW FILL Mode and then into the Stasis-PURGE Mode, followed by POST-PURGE FILL Mode.



NOTE

After Post-PURGE FILL Mode, if any additional STASIS Cycles are programmed, the 931 will re-enter STASIS Mode.


8. Upon completion of the last preset STASIS Mode cycle, the 931 will then advance into the COOL Mode, SLOW FILL Mode, PURGE Mode and POST-PURGE FILL Modes sequentially just as if it were in an Auto Mode.



NOTE:

Chamber Lid may not be secured properly or damaged O-Ring could cause chamber pressure to drop during STASIS Mode.


- During STASIS Mode, if chamber pressure drops to 0 PSI, the 931 will pause and display a Chamber Security Check Sign (Security Feature) stating “Confirm Chamber is Secured”.




NOTE:

Please refer to “Auto Mode” (Chapter 4.2) for detailed instructions.

- At this point until the end of the process run, the 931 will process each Mode Automatically.



POST-PURGE Mode




HEAT Mode




NOTE


tousimis-Equilibrium Time (CP Time) may need to be extended as per Thermal Conductivity of the sample.



tousimis Equilibrium



BLEED Mode




VENT Mode



NOTE:
If LCO₂ Tank runs out during a process run, press **PAUSE** and replace the LCO₂ Tank.
Press **RESUME** to reactivate process run after the new LCO₂ Tank has been installed.


4.5 Pause & Resume During Process Run

1. At any time during the Auto, Manual, or STASIS Modes the **PAUSE** button will allow the 931 to enter into the **PAUSE** Mode.



The 931 can then resume the process by pressing the **RESUME** button.


2. The **MODIFY PROCESS SETTINGS** button will redirect the display view to the “Modify Settings” page where the 931 Auto Mode Settings may be modified.




See instructions in “Auto Mode” (Chapter 4.2) to modify Auto Mode Settings and press **RESUME** to continue process run.

Pressing **MODIFY PROCESS SETTINGS** during “Manual Mode” will advance 931 from “Manual Mode” to “Auto Mode.” The remaining process run steps will **RESUME** as Auto Mode.

3. During STASIS Mode, there will be an extra button for STASIS Setting Modification. See “STASIS Mode” (Chapter 4.4) instructions to modify settings.



4. Press **QUIT PROCESS** button on the Pause Screen and Modify Settings Screen will display a **BLEED PORTS** and **VENT PORTS** option to reduce chamber pressure in order to terminate the current process run.



BLEED PORTS: This vents the chamber pressure back to 0 PSI while allowing intermediary fluid (Ultrapure Alcohol) to stay within the Process Chamber. Therefore, your sample will not be exposed to gas.

VENT PORTS: This vents the chamber pressure out and it is often set at a quicker flow. However, all intermediary fluid (Ultrapure Alcohol, if present) in the Process Chamber will be discharged through these ports.



4.6 Recipe Mode

1. The “Recipe Mode” allows the input and storage of up to 10 custom process run protocols. Press **RECIPE** button from “Main Menu” page to enter Recipe Mode. To create a new recipe, press **ENTER NEW RECIPE** button. To run or modify a previously saved recipe, press previously saved recipe name.

2. Press **EDIT** button next to either “Recipe Name”, “Recipe Owner” or “Sample Name” to update the appropriate information for the Recipe.

* Press **Stasis Mode** Toggle Switch to set **Stasis Mode** ON/OFF for the recipe.

3. Type "Recipe Name" and press **Enter** to update information.

4. Press **NEXT** when the complete information has been entered.
Press **MAIN MENU** or **RECIPE** button to exit without saving.

5. If the **Stasis Mode Switch** is **OFF**:

You may modify the preset values. Pressing **START** will save recipe and start the "Auto Mode" immediately. **SAVE RECIPE** will Save the recipe only.

Recipe Name	Hanz Tp 334 r1	
START	01 : 30	Slow FILL
SAVE RECIPE	02 : 30	FILL Time
RECIPE	10 : 00	PURGE Time
MAIN MENU	01 : 30	Post PURGE
	02 : 00	Critical Point
	550 PSI	BLEED/VENT

Press **RECIPE** and **MAIN MENU** to exit without saving.

6. If the Stasis Mode Switch is ON:

You may modify the preset values and press **NEXT**. Press **RECIPE** and **MAIN MENU** to exit without saving.



The **Recipe Setup** screen displays the recipe name **Hanz Tp 334 r1**. On the left, there are three navigation buttons: **NEXT**, **RECIPE**, and **MAIN MENU**. On the right, a list of parameters and their values is shown, each with a corresponding button: **01 : 30** (Slow FILL), **02 : 30** (FILL Time), **10 : 00** (PURGE Time), **01 : 30** (Post PURGE), **02 : 00** (Critical Point), and **550 PSI** (BLEED/VENT).

7. Set “Stasis Time” and “Stasis Cycle” (See Chapter 4.4 “Stasis Mode”) and press **START** to save recipe and start “Stasis Mode” immediately.

Press **SAVE** to save the recipe only.




The **Recipe Stasis Setup** screen displays the recipe name **Hanz Tp 334 r1**. Below the name, it shows the current settings: **Stasis Time (HH:MM)** as **02 : 00**, **Stasis Cycles** as **03**, and **Stasis PURGE** as **03**. Each value is displayed above a set of up/down arrow buttons for adjustment. Below these are buttons for **Hours**, **Minutes**, **Cycles**, and **Cycles**. On the right side, there are three buttons: **START**, **SAVE**, and **UNDO**. At the bottom, there are two navigation buttons: **RECIPE** and **GO BACK**.

Press **GO BACK** and **RECIPE** to exit without saving.


4.7 Critical Point Failure Correction

Symptom: The “CP Failure” sign will be displayed if the Temperature and/or Pressure in the Process Chamber has not reached the “Critical Point” in the proper sequence during HEAT Mode.



A successful “Critical Point” process (tousimis Equilibrium) is achieved when the internal micro processor receives the following two signals in this sequence:

1. Pressure Signal sent via High Pressure Sensor between 1150-1200 PSI.
2. Heat Signal sent via Heat Thermostat shutting OFF typically between 34-38°C



NOTE:

- + : CP Temperature
- : CP Pressure
- ⊕ : Critical Point

The + sign on the Pressure/Temperature display indicates that the Process Chamber has reached critical point temperature. An ○ sign indicates that the Process Chamber has reached critical point pressure.


Critical Point Failure Correction Options:

Provided the micro processor receives the signals in the previously mentioned sequence, there will not be a 'Critical Point Failure.'


However, should the micro processor receive these signals in the reverse order - Heat Signal (1st) followed by Pressure Signal (2nd) - a failure will be indicated by 'CP Failure' flashing.

Solutions:

Scenario #1: Pressure is greater than 1072 PSI.




Press **BLEED**. 'CP Failure' display will be turned off and the 931 will default back to the normal Auto-Mode and resume in **BLEED** Mode moving on to completion as normal.



The LCO₂ tank should be replaced with a new LCO₂ tank prior to commencing the next process run.

Scenario #1 may be deployed as the theoretical physical completion of the critical point (1072 PSI + 31°C) for CO₂ has been attained even though the system's sensors did not detect it. The system's sensors are designed to over-ride both pressure and temperature minimums for the CO₂ critical point.


Scenario #2: Pressure is less than 1072 PSI.



Replace LCO₂ tank with a new LCO₂ tank.

DO NOT OPEN CHAMBER AT THIS POINT.
KEEP SAMPLE PRODUCT WITHIN PROCESS CHAMBER.

Press **COOL** button to re-initiate the process.




Press **SLOW FILL** button as per normal initiation of a new process run.



4.8 Power Failure Recovery

Symptom: Power outage during process run.

1. Power ON when the facility power is back ON.
Leave VENT Valve switch in CLOSE position.



Press **Enter User/Sample Information** button to enter User/Sample Information.

Press **Main Menu** button to start operation without setup User/Sample Information.

2. Press **MANUAL** button to initiate "Manual Mode" process.



3. Resume process by pressing the button of the process mode (COOL, SLOW FILL, PURGE, HEAT, BLEED, VENT) at which the process had stopped during power outage.







NEW USER

4.9 New User

1. To update User & Sample Information, press **NEW USER** button.



2. Enter User Name or ID (maximum 18 characters) and press **Enter**.



3. In the following screen, Enter Sample ID or Description (maximum 18 characters) and press **Enter**.

Please Enter Sample Description

—


1	2	3	4	5	6	7	8	9	0	Clr
Q	W	E	R	T	Y	U	I	O	P	
A	S	D	F	G	H	J	K	L		
Z	X	C	V	B	N	M	Backspace			
SKIP		Space						Enter		




SETTINGS

4.10 Settings (Date/Time)


1. To update Date/Time information, press **DATE/TIME** button.



2. Select Calendar Display Type.




3. Select Clock Type.



The screen displays the title "Clock Type" in a dark header bar. To the right of the title, the text "Hanz Wetzler", "Tp 334 Btch 16653", and "01/26/2012 10:43 AM" is visible. Below the header, the current pressure "0 PSI" and temperature "15°C" are shown. The main instruction "Select Clock Type" is centered. Two buttons, "12 Hrs" and "24 Hrs", are presented for selection. A "GO BACK" button with a left-pointing arrow is located at the bottom left.


4. Using arrow buttons below each time variable, adjust time. Press **Save** to update time. **Reset** to reset the time. Press **Exit** to go back to "Settings" menu.




The screen displays the title "Set Date & Time" in a dark header bar. To the right of the title, the text "Hanz Wetzler", "Tp 334 Btch 16653", and "01/26/2012 10:43 AM" is visible. Below the header, the date and time are displayed in a grid: "Month: 02", "Date: 27", "Year: 2012" in the top row, and "Hour: 09", "Minute: 23", and "AM" in the bottom row. Each value is flanked by left and right arrow buttons for adjustment. Below the date and time fields, there are three buttons: "Reset", "Exit", and "Save".

4.11 Flow Charts

Auto Mode



Stasis Mode

Recipe Mode

